

source



UG001 SCIIN One

User Instructions

Order yours today:

T. 01276 676 006

E. info@sourcesupplies.co.uk

W. sourcesupplies.co.uk

reduce + reuse + recycle

SCIIN ONE

Antiviral Face Mask

UG001 / 01/07/2020

User instructions

How to wear

Hold the mask with the curved edge at the bottom and simply hook over your ears (see diagram).

Ensure your nose, mouth and chin are covered for maximum protection.

Care guidelines

We recommend hand-washing your mask for up to 15 washes.

Do not tumble dry as this will cause the material to shrink and damage the antiviral coating. Allow to air dry naturally on a flat surface.

To remove creases you can iron the product on a low settings (silk) making sure the logo is face down and no steam is used.

Keep your mask in a cool dry place when not in use and/or in the original packaging.

Best practise

Avoid touching your face after removing your mask and always wash your hands to avoid transferring any germs.

Disposal




The material is 100% biodegradable and safe to be disposed of in a normal waste bin. The packaging is 100% recyclable. Please dispose of responsibly to protect the environment.

Antiviral coating

The antiviral coating provides 99% effective protection from new. After 15 hand-washes the coating is still 80%+ effective according to Polygiene's tests. Polygiene's antimicrobial technology prevents growth of microorganisms in the fibres (such as bacteria, fungi and parasites) maintaining a cleaner fabric for longer, and the non-necessity for frequent washing.

Allergies

ViralOff is safe and has no detrimental effect on your skin nor does it interfere with the skin's natural bacterial flora.



 **ViralOff**
POLYGIENE TECHNOLOGY

Order yours today: T. 01276 676 006 E. info@sourcesupplies.co.uk W. sourcesupplies.co.uk

SOURCE

Frequently Asked Questions

Can I wear the Sciin One both ways?

Yes. The Sciin One has coating on both sides so it doesn't reduce the antiviral effect whichever way you wear it.

Does the Sciin One lose its shape?

The Sciin One is made from lightweight bamboo viscose for increased comfort and breathability so you will see a slight stretch over time depending on how long you wear the mask. Hand washing and allowing to dry naturally on a flat surface will help return the mask back to its original shape and size.

Can a child wear a Sciin One?

Sciin One is available in child sizes on request, please email us at hello@sciin.co.uk to enquire.

Can I wash my Sciin One?

We recommend a maximum of 15 hand washes in warm water without detergent.

How often should I wash my Sciin One?

You can wear your Sciin One for as long as you like before having to wash it and you do not have to wash the Sciin One after each wear, however we do recommend washing if the fabric becomes dirty. The antiviral coating provides 99% effective protection from new. After 15 hand-washes the coating is still 80%+ effective.

Can I dry my Sciin One on a radiator or in a tumble dryer?

We do not advise drying on a radiator or in a tumble dryer and recommend leaving to air dry naturally on a flat surface. High temperatures can damage the antiviral coating and reduce its effectiveness.

Can I iron my Sciin One?

We do not recommend ironing Sciin One but if you wish to remove any creases we recommend ironing with the logo face down on the lowest setting (silk) and no steam is used. High temperatures can damage its antiviral coating.

Can I dry clean my Sciin One?

No. We do not recommend dry cleaning your Sciin One.

Will it provide 100% protection against viruses?

No device can guarantee 100% protection against viruses. The antiviral coating provides 99% effective protection from new. After 15 hand-washes the coating is still 80%+ effective.

How else can I protect myself against viruses?

Wash your hands frequently with soap and clean water. Cough and sneeze into a tissue and dispose of it immediately. Avoid touching your eyes, ears, nose and mouth.