

URIZAP Shock

Technical Data Sheet



Intensive uric scale removal for problematic urinals

Unique oxygen-based formulation dissolves uric acid and slime

Sources of odour and causes of blockages are rapidly eliminated

Pair with **URIZAP** for complete urinal maintenance solution

USAGE INSTRUCTIONS:

- 1. Pour 1 (20 – 25 g) scoop of URIZAP Shock directly into the urinal drain**
 - Rinse down with **half a cup (~ 150 mL) of water**. Warm water works best
 - The **effervescent action** activates the formulation. Depending on urinal plumbing design, **a small quantity of foam may appear in the drain.**
 - Leave **URIZAP Shock** in place for **at least 1 hour, preferably overnight.**
 - Formulation is pH-neutral, any residue can be safely rinsed away
 - Repeat dosing **daily for 1 week** to eliminate uric scale, slime and buildup
 - Follow up with regular **URIZAP** dosing for deodorization and maintenance
-

PRODUCT INFORMATION:

Application areas:	Urinals – ceramic, stainless steel, resin	
Product size:	1000 g	40 x 25 g doses per pack
Shelf life:	12 months	From date of manufacture
Product appearance:	Blue	Granulated powder
Effective temp. range:	10 – 50 °C	
Product pH:	7 – 9	Neutral, non-corrosive
Product fragrance:	Menthol	
Active component:	Peracetic acid	Organic oxidant, dissolves uric salts, biofilm and slime
Hazard class (Granules):	H318	Causes serious eye damage
Recommended PPE:	Standard equipment for housekeeping services Wash hands after use	